
# IKA

designed for scientists

## RCT basic *safety control*

Betriebsanleitung Ursprungssprache	DE	<b>6</b>
Operating instructions	EN	<b>17</b>
Mode d'emploi	FR	<b>28</b>
Руководство по эксплуатации	RU	<b>39</b>
Instrucciones de manejo	ES	<b>50</b>
Instruções de serviço	PT	<b>61</b>
Návod na použitie	SK	<b>72</b>
使用说明	ZH	<b>83</b>
取扱説明書	JA	<b>94</b>
사용 설명서	KO	<b>105</b>
Veiligheidsaanwijzingen	NL	<b>116</b>
Avvertenze per la sicurezza	IT	<b>118</b>
Säkerhetsanvisningar	SV	<b>120</b>
Sikkerhedshenvisninger	DA	<b>122</b>
Sikkerhetsveiledning	NO	<b>124</b>
Turvallisuusohjeet	FI	<b>126</b>
Wskazówki bezpieczeństwa	PL	<b>128</b>
Bezpečnostní pokyny	CS	<b>130</b>
Biztonsági utasítások	HU	<b>132</b>
Varnostna navodila	SL	<b>134</b>
Ohutusjuhised	ET	<b>136</b>
Drošības norādījumi	LV	<b>138</b>
Saugos nurodymai	LT	<b>140</b>
Указания за безопасност	BG	<b>142</b>
Indicații de siguranță	RO	<b>144</b>
Υποδείξεις ασφάλειας	EL	<b>146</b>

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## Contents

	Page
EU Declaration of conformity	17
Safety instructions	17
Intended use	19
Unpacking	19
Assembling the stand	19
Operation	20
Setting operation mode	21
Er5	21
Setting the safe temperature limit	22
Controlling the medium temperature limit using a contact thermometer	23
Setting the temperature control mode	24
Maintenance and cleaning	24
Accessories	24
Error codes	25
Technical Data	26
Warranty	27

## EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at [sales@ika.com](mailto:sales@ika.com).

### /// General information

- **Read the operating instructions in its entirety before using the device and follow the safety instructions.**
- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Electrical outlet must be grounded (protective ground contact).
- **Attention – Magnetism!**  
Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).
- **Risk of burns!**  
Exercise caution when touching parts of the housing and the heating plate.  
The heating plate can reach dangerous temperatures. Pay attention to the residual heat on the heating plate after switching off the stirrer. The device may only be transported when the heating plate has cooled down.

### /// Device design

- Do not use the device in explosive atmospheres, it is not EX-protected.
- With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the device must be clean and undamaged.
- Ensure that the power cord set / temperature sensor cable does not touch the heating plate.
- Check the device and accessories for damage before each use. Do not use damaged components.

### /// Permissible medium / contaminants / side reactions

- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any ex-

tra energy produced in other ways, e.g. through light irradiation.

- Beware of hazards due to:
  - flammable materials,
  - combustible media with a low boiling temperature,
  - glass breakage,
  - incorrect container size,
  - overfilling of media,
  - unsafe condition of container.
- Process pathogenic materials only in closed vessels under a suitable fume hood.
- In unsupervised and safe operation, this device may only handle or heat media whose flash point is above the set safety temperature limit. The safety temperature limit must always be set to at least 25 °C lower than the flash point of the media used. (acc. to EN 61010-2-010)
- The heating plate can heat up due to the action of the magnets at high motor speeds, even if the heater is switched off.
- Please consider any possible contaminations and unwanted chemical reactions.
- It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C – 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.  
(Source: *Römpps Chemie-Lexikon* and *"Ulmann"*, Volume 19)

#### /// Procedures during sample runs

- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. There may be a risk from:
  - splashing and evaporation of liquids,
  - ejection of parts,
  - release of toxic or combustible gases.
- Reduce speed if:
  - medium splashes out of vessel because the speed is too high,
  - device is not running smoothly,
  - container moves on the base plate,
  - an error message is displayed.

#### /// Accessories


- Safe operation is guaranteed only with the use of original IKA accessories.
- Ensure that the external temperature sensor is inserted into the medium to a depth of at least 20 mm when connected.
- Always disconnect the plug before attaching accessories.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Observe the operating instructions of the accessories.

#### /// Power supply / Switching off the device

- The device will automatically restart in mode b following any interruption to the power supply.
- The voltage stated on the type plate must correspond to the mains voltage.
- The outlet for the mains plug must be easily accessible.
- The device can only be disconnected from the mains outlet by pulling out the power cord set or the connector plug.

#### /// For protection of the equipment

- The device may only be opened by qualified and IKA approved experts.
- Do not cover the device, even partially e.g. with metallic plates or film. This may result in overheating.
- Protect the device and accessories from bumps and impacts.
- Keep the base plate clean.
- Observe the minimum distances:
  - between devices min. 100 mm,
  - between device and wall min. 100 mm,
  - above the device min. 800 mm.

## Intended use

/// Use

- The magnetic stirrer is suitable for mixing and / or heating substances.

/// Area of use

- Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- The safety of the user cannot be guaranteed:
  - if the device is operated with accessories that are not supplied or recommended by the manufacturer,
  - if the device is operated improperly or contrary to the manufacturer's specifications,
  - if the device or the printed circuit board are modified by third parties.

## Unpacking

/// Unpacking


- Please unpack the device carefully
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder)

/// Scope of delivery

- RCT basic
- Mains cable
- Operating instructions
- Screwdriver
- Protective cover
- Temperature sensor PT 1000

## Assembling the stand

- › Review the mounting and safety instructions of the IKA boss head clamp prior to using it. (see "Accessories")
- › The device must not be suspended from the support rod!


**Note:** When using bath attachments with a diameter over 180 mm, use support rod in conjunction with an extension (see "Accessories").

### Risk of tipping!

Make sure that the center of mass of the attached device does not protrude over the safe area indicated below with a dotted rectangle.

## Operation

Ensure that the protective film is removed from the base plate before use!

<b>Commissioning</b>	<ul style="list-style-type: none"> <li>☞ Move mains switch (A) to the OFF position.</li> <li>☞ Insert the mains power cable into the power socket (K).</li> <li>☞ Move mains switch (A) to the ON position.</li> <li>⇒ The unit will be set to factory setting mode A (see "Setting operation mode").</li> </ul>					
<b>Stirring</b>	<ul style="list-style-type: none"> <li>☞ Adjust the motor speed using the rotating / pressing knob (D).</li> <li>⇒ The value selected will be shown on the display (F).</li> <li>☞ Set the start point for the agitation function by pressing the rotating / pressing knob (D).</li> </ul> <p>① ⇒ The displayed value will flash until the desired motor speed is reached.</p>					
<b>Heating</b>	<ul style="list-style-type: none"> <li>☞ Set the safe temperature limit (see "Setting the safe temperature limit").</li> <li>☞ Adjust the set-point temperature using the rotating / pressing knob (C).</li> <li>⇒ The value selected will be shown on the display (E).</li> <li>☞ Set the temperature control mode (see "Setting the temperature control mode").</li> <li>☞ Set the start point for the heating function by pressing the rotating / pressing knob (C), also on use of a contact thermometer.</li> </ul> <p>The set point and actual temperatures will be shown alternately on the display (E).</p> <p>① ⇒ When the heating is switched on, the LED "Heating plate" (G) and the LED "Set point value" (I) will be lit.</p> <p style="border: 1px solid black; padding: 2px; display: inline-block;">The maximum temperature that can be set for the heating plate is 310 °C.</p> <p>① ⇒ During agitation and standby operation, the display (E) will show <b>H0E</b> if the plate temperature is above 50 °C.</p>					
<b>Connecting external thermometers</b> <i>(direct temperature control in the media)</i>	<ul style="list-style-type: none"> <li>☞ Move mains switch (A) to the OFF position.</li> <li>☞ Detach contact plug (L).</li> <li>☞ Attach a DIN 12878 (Class 2) compliant safety contact thermometer or a PT 1000 temperature sensor to the connection (L).</li> <li>☞ Move mains switch (A) to the ON position.</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <tbody> <tr> <td style="width: 35%; padding: 5px;">① Temperature sensor PT 1000</td> <td style="padding: 5px;">⇒ The actual temperature for the temperature sensor shown on display (E) will correspond to the temperature of the media. "LED, external temperature sensor" (H) will be lit.</td> </tr> <tr> <td style="padding: 5px;">① Contact thermometer e.g. ETS-D5</td> <td style="padding: 5px;">⇒ Follow the operating instructions for the contact thermometer. "LED, external temperature sensor (decimal point for the temperature display)" (J) will be lit. ⇒ When a contact thermometer is connected, the display (E) will only show the set-point temperature that has been set.</td> </tr> </tbody> </table>		① Temperature sensor PT 1000	⇒ The actual temperature for the temperature sensor shown on display (E) will correspond to the temperature of the media. "LED, external temperature sensor" (H) will be lit.	① Contact thermometer e.g. ETS-D5	⇒ Follow the operating instructions for the contact thermometer. "LED, external temperature sensor (decimal point for the temperature display)" (J) will be lit. ⇒ When a contact thermometer is connected, the display (E) will only show the set-point temperature that has been set.
① Temperature sensor PT 1000	⇒ The actual temperature for the temperature sensor shown on display (E) will correspond to the temperature of the media. "LED, external temperature sensor" (H) will be lit.					
① Contact thermometer e.g. ETS-D5	⇒ Follow the operating instructions for the contact thermometer. "LED, external temperature sensor (decimal point for the temperature display)" (J) will be lit. ⇒ When a contact thermometer is connected, the display (E) will only show the set-point temperature that has been set.					

## Setting operation mode

### Operating the device in mode A, b or d

#### Mode A

All settings will be stored if the device is switched off or disconnected from the power supply. The agitation and heating functions will be set to OFF when the device is powered on.

#### Mode b

All settings will be stored if the device is switched off or disconnected from the power supply. The agitation and heating functions will be set to ON or OFF when the device is powered on, depending on the previous status of the device.

#### Mode d

In Mode d the device behaves the same as in Mode A with the exception that:

- The temperature setting must be confirmed by pressing the temperature knob. To change the temperature setting, turn the temperature knob until the desired value appears.
- The new value flashes for 5 seconds on the display. Confirm the new temperature setting by pressing the temperature knob, otherwise the temperature setting value will skip back to its previous value.

### *Factory setting: mode A*

The mode selected will be shown on the display when the device is started up.

### Change the mode

- ☞ Move mains switch (A) to the OFF position
- ☞ Press and hold rotating / pressing knobs (C and D)
- ☞ Move mains switch (A) to the ON position
- ☞ Release rotating / pressing knobs (C and D)
  - ⇒ The set value is indicated on the display (F)

xx A

**Note:** Mode d is available from software version 40

## Er5

Providing safety protection, Er5 is taken to recognise that the sensor has not been placed in the medium when the heating is switched on.

**Note:** The user can set a value from 1 to 30 min for this time limit depending on the application.


If time limit is set to 0 min, the Er5 detection is disabled.

This function will only be active if:

- Sensor temperature is  $< 50\text{ }^{\circ}\text{C}$
- Difference target temperature / sensor temperature  $> 5\text{ K}$

### *Factory setting: 5 min*

The time selected will be shown on the display when the device is started up if the sensor is connected.

### Changing the Er5 settings

- ☞ Set target temperature to  $5\text{ }^{\circ}\text{C}$
- ☞ Switch off device
- ☞ Hold down temp knob
- ☞ Switch on device
- ☞ Select time value from 0...30 min (in 1 min steps) by turning the speed knob
- ☞ Release the temp knob

**Note:** Er5 is available from software version 40.

## Setting the safe temperature limit

The maximum achievable heating plate temperature is restricted by an adjustable safety temperature limit. Once this limit has been reached, the device stops heating.

**The safe temperature limit must always be set at least 25 °C lower than the fire point of the media to be processed! (acc. to EN 61010-2-010)**

The adjustable maximum heating plate temperature is at least 10 °C under the set safety temperature limit.

Setting range: [50 °C] to [max. SET temperature + 50 °C]

Factory setting: [max. SET temperature + 50 °C]

### Setting the safety temperature limit

After switching on the device, the safety temperature limit can be adjusted using a screwdriver.

Do not turn the setting screw beyond the clockwise or anticlockwise stop. This will cause irreparable damage to the potentiometer.

- Set the main switch to the ON position.
- Using the screwdriver supplied, turn the “Safe Temp” setting screw to the clockwise stop.
- Use the temperature rotating / pressing knob to set the target temperature to the desired “Safe Temperature” and wait until the temperature is achieved.
- Turn the “Safe Temperature” setting screw slowly anticlockwise until the heating function switches off and the display shows **Er25**.
- The “Safe Temp” value is displayed on the screen.

Setting the safety temperature limit has to be done as described. The value “Safe Temperature” shown on the display only serves for visualization.

### Functional test: safety circuit switch-off

- Heat the device to a temperature of over 50 °C.
- Set the safety temperature limit as far to the left as possible (50 °C) and set the mains switch to the ON position.
- The display will show: **Er25**

## Controlling the medium temperature limit using a contact thermometer

The preferable method for controlling the average temperature is with contact thermometer. After the set point temperature has been adjusted, this results in a short heating-up time, practically no temperature drift and only minor fluctuation in temperature.

A 6-pin jack is located on the rear side of the device for connecting the PT 1000 series, contact thermometer or the contact plug. The electronics of the devices returns a test current that must flow via connector pins 3 and 5 for the heating plate to heat up.

### **Safety contact thermometers:**

Acc. to DIN 12 878 class 2 or acc. Gerstel are connected with a three-wire cable, the test current flows through the contact thermometer.

### **Safety function:**

If the test current is interrupted because of e.g. breakage of contact thermometer or falling out of the cable plug, the heating cuts off.

### **Contact thermometer without safety circuit:**

Acc. to DIN 12 878 class 0. The device only heats if the test current circuit is closed by an electrical connection of the plug pins 3 and 5.

### **2-wire connecting cables:**

Connect plug pins 3 and 5 of the device plug.

### **3-wire connecting cables:**

Here the test current circuit can also be produced in the terminal head of the contact thermometer (connect plug pins 2 and 3). A 3-wire cable with the required bridge is available (accessories).

### **Settings:**

For detailed instructions for settings and limit values, please refer to the operating manual of the device you are connecting.

**The desired average temperature can be adjusted on the contact thermometer. The required surface temperature of the heating plate can be selected with the rotating / pressing knob or button.**

Adjusting the temperature of device to the maximum adjustable temperature will result in the fastest possible heating time, but the average temperature may fluctuate to values above the set-point temperature on the contact thermometer. By adjusting the rotating / pressing knob or button to approximately twice the set-point value of contact thermometer (with a set-point of + 60 °C, the temperature of device would be set to + 120 °C), you will reach a good compromise between a fast heating time and over-shooting the set point. If you adjust the temperature of the device to exactly the set-point temperature, the medium will not reach the set-point temperature because some loss of the heat will always occur between the heating plate and the medium.

**If a fault occurs in the control circuit, the maximum heating plate temperature can be adjusted with the safety temperature setting screw.**

## Setting the temperature control mode

When using an external PT 1000 temperature sensor, the user can choose between two types of control:

### PI mode

Good control results, minimized overshooting, slow rise in temperature

### 2P-mode (two-point controller)

Maximum heating rate, increased overshooting

*Factory setting: PI mode*

### Changing the temperature control mode

1. Set target temperature to 2 °C
2. Switch off device
3. Hold down temp knob
4. Switch on device

This activates or deactivates the two-point controller (2P) for the external PT 1000 temperature sensor (2P appears on the display).

When switching on the device, 2P is also displayed next to the operating mode A/b/d when the two-point controller is activated. When the heating function is active in 2P mode, (with external PT 1000 temperature sensor) the actual temperature flashes - that notifies the user that the process temperature has been greatly overshoot!

- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

### Ordering spare parts

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item number and designation of the spare part, see **www.ika.com**.
- Software version

### Repair

**Please send instrument in for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.**

For this you should request the “**Safety Declaration (Decontamination Certificate)**” from **IKA**, or use the download printout of it from the **IKA** website **www.ika.com**.

Return the instrument in its original packaging. Storage packaging is not sufficient. Also, please use suitable shipping package materials.

## Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

### Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean IKA devices.  
These are water (with tenside) and isopropanol.

## Accessories

- For further accessories see **www.ika.com**.

## Error codes

Any malfunctions during operation will be identified by an error message on the display (E and F).  
 Proceed as follows in such cases:

- Switch off device using the main switch (A)
- Carry out corrective measures
- Restart device

Error code	Cause	Effect	Solution
Er3	Temperature inside device too high	Heating off	- Switch off device and allow to cool down
Er4	Motor blockage	Heating off Motor off	- Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the motor</i>
Er5	No temperature increase in sensor at continuous temperature difference	Heating off	- Place the sensor in the medium - Reduce the volume of the media - Use a carrier fluid with better heat conductivity properties - Replace the glass vessel with a metal pot - Increase the "Time-out" period
Er6	Break in safety circuit	Heating off	- Plug in contact plug (L) - Plug in PT 1000 contact thermometer / temperature sensor - Replace faulty connecting cable, plug, or contact thermometer
Er24	Surface temperature ( <b>temperature of control sensor</b> ): of the base plate is higher than the set safe temperature limit	Heating off	- Switch off device until the surface temperature of the base plate is lower than the selected safe temperature limit - Set a higher safe temperature limit
Er44	Surface temperature ( <b>temperature of safety sensor</b> ): of the base plate is higher than the set safe temperature limit	Heating off	- Switch off device until the surface temperature of the base plate is lower than the selected safe temperature limit - Set a higher safe temperature limit
Er25	Heating and switching element monitoring	Heating off	- Switch off device - Set the safe temperature limit > 55 °C - see also "Functional check of inactivating the safety circuit" - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the heating element</i>
Er26	Difference between temperature of safety sensor and temperature of control sensor <b>control temperature</b> > (safety temperature + 40 K)	Heating off	- Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i>
Er46	Difference between temperature of safety sensor and temperature of control sensor <b>safety temperature</b> > (control temperature + 40 K)	Heating off	- Switch off device - <i>Warning! Only to be carried out by authorized service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i>

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault.

## Technical Data

### Device

Operating voltage range	Vac	(220 ... 230) ± 10%
	Vac	115 ± 10%
	Vac	100 ± 10%
Nominal voltage	Vac	230/ 50 Hz
	Vac	115/ 60 Hz
	Vac	100/ 60 Hz
Frequency	Hz	50/ 60
Power consumption (+10%) max. at	230 Vac	650
	115 Vac	610
	100 Vac	610
Display		LED
Permissible duration of operation	%	100
Permissible ambient temperature	°C	+ 5 ... + 40
Permissible relative humidity	%	80
Protection type acc. DIN EN 60529		IP 42
Protection class		I
Overvoltage category		II
Contamination level		2
Operation at a terrestrial altitude	m	max. 2000
Dimensions (W x D x H)	mm	160 x 270 x 85
Weight	kg	2.5
<b>Motor</b>		
Number of stirring positions		1
Direction of rotation		right
Speed control		rotating / pressing knob
Speed range	rpm	50 ... 1500
Motor rating output	W	9
Setting resolution	rpm	10
Stirring bar length	mm	20 ... 80

Speed variation (no load, nominal voltage, at 1500 rpm, ambient temperature 25 °C)	%	± 2
Self-heating of the heating plate by max. stirring (RT:22 °C / duration:1 h)	K	+ 28
Stirred quantity max. (H <sub>2</sub> O)	ltr	20
<b>Base plate</b>		
Dimensions	mm	Ø135
Material		Al-alloy
<b>Heating</b>		
Heating power (-5%/+10%) at nominal voltage	W	600
Adjustment and display resolution	K	1
Surface temperature	°C	Ambient temperature ... 310
Heat control		rotating / pressing knob
Heating rate heating plate	K/min	6.5
Temperature sensor PT 1000 variation DIN EN 60751 Kl. A	K	≤ ± (0.15 + 0.002x T )
Temperature variation max. at 100 °C	K	± 1.5
Heating plate control hysteresis, no container, center of heating plate at 100 °C	K	± 5
Control hysteresis with temperature sensor PT 1000*	K	± 1
with ETS-D5*	°C	± 0.5
with ETS-D6*	°C	± 0.2
Sensor in medium detection (Error 5)		yes
<b>Adjustable safety circuit</b>		
Safety temperature limit (adjustable)	°C	(50 ... 360) ± 10

\* The control accuracy values given were determined using the following set-up: 500 ml water in 600 ml glass beaker, magnetic stirring bar 40, 600 rpm, 50 °C  
Subject to technical changes!

## Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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