Enthalpy and entropy as driving forces of chemical reactions

(Item No.: P1136300)

Curricular Relevance

Difficulty

Difficult

Preparation Time

10 Minutes

Execution Time

20 Minutes

Recommended Group Size

1 Student

Additional Requirements:

Experiment Variations:

Keywords:

Principle and equipment

Principle

The first experiment described here can be used to introduce the concept of enthalpy as a driving force of chemical reactions. High-energy chemical systems tend to assume a lower-energy, more stable state or to enter a lower-energy system while releasing energy. This becomes particularly clear in the case of reactions that occur and proceed spontaneously without the supply of activation energy (ignition energy, start energy). The second experiment can be used to introduce the concept of entropy as a driving force of chemical reactions. The transition from a solid, ordered state to a less ordered state (liquid) or to a completely disordered state is usually accompanied by the absorption of energy (see also the phenomenon of enthalpy of solution with a reduction in temperature (CT 24.5) during the dissolution of many solid substances in solvents or the absorption of heat during melting and evaporation processes). Evaporation is a good example of the tendency of substances to transition to the gaseous state even at normal temperatures

The aim of the experiments in chapter 11 (energy turnover during chemical reactions) is the introduction of reaction processes. Stable and metastable systems are used to explain the effect of activation energy and catalysts. Once the release of energy has been recognised as an important driving force of the reactions, the enthalpies can be determined in a quantitative manner in order to understand the reaction processes. An important result is the realisation that the energy turnover of a reaction can be measured relatively easily.

Tasks

1. Spontaneous exothermic reactions
2. Spontaneous endothermic reactions

Equipment

Safety information

Concentrated acids are highly caustic. They burn the skin and destroy textile fabrics. For diluting them, first add the water and then the acid (protective glasses, laboratory coat, gloves).

First aid

Rinse the affected skin areas and eyes with the lid gap wide open thoroughly with plenty of water.

Disposal

Solutions must be diluted with water, neutralised (pH 6-8), and flushed away. Collect any solutions that contain heavy-metal ions in a collecting vessel for heavy-metal salt solutions.

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
Observation

Task 1

Strong formation of hydrogen and a temperature rise of the solution by approximately 10°C.

Reaction equation:

After this observation, the experiment is repeated with an equivalent mass of magnesium (0.75 g).

Temperature rise by several degrees.

Task 2

The mixture of solid substances quickly turns into a solution (liquid). At the same time, the smell of ammonia can be noticed. After approximately 2 minutes, the beaker has frozen to the supporting block. The temperature measuring device indicates a temperature of several degrees below zero (up to -20°C).

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