


The effect of acids and bases on natural and synthetic indicators

(Item No.: P7510200)

Curricular Relevance



Difficulty



Easy

Preparation Time



10 Minutes

Execution Time



20 Minutes

Recommended Group Size



2 Students

Additional Requirements:

- Beetroot dye
- Flower dye
- Natural dye (red cabbage juice)

Experiment Variations:

Keywords:

acid, base, indicator, pH-value

Information for teachers

Introduction

Application

Synthetic indicators can aid in the identification of unknown substances, for example first classifications of acids and bases.



Experiment set-up

Educational objectives

Indicators are dyes, which identify acids or bases with a change of colour. Indicators can be either synthesized or natural dyes.

Task

The students will gain first experience with indicators and how they work.

Prior knowledge

This experiment serves as an introduction into the topic 'Acids, bases and indicators'. No prior knowledge required.

Principle

Indicators are organic dyes. They show structural changes with a varying pH value. The structural change in the molecule is responsible for the change in colour.

Notes concerning the experiment and set-up

Prepare the following solutions:

5% hydrochloric acid (Add 50 ml distilled water to a suitable volumetric flask, pipette 13 ml of hydrochloric acid, 37 % and fill up to 100 ml with distilled water)

10% acetic acid (Add 50 ml distilled water to a suitable volumetric flask, pipette 11 ml of concentrated acetic acid and fill up to 100 ml with distilled water)

5% NaOH (Dissolve 5 g NaOH in 95 g distilled water)

5% KOH (Dissolve 5 g KOH in 95 g distilled water)

Prepare the following natural indicators (dyes):

Beetroot dye: Cut a beetroot into pieces and cover them with hot water. Wait 10 minutes and filter the solution.

Red Cabbage/flower: Add 4 spoons of sand and 10 ml of methylated spirit with red cabbage/flower to a mortar. Grind with a pestle. Add another 5 ml of methylated spirit. Grind for 5 minutes. Add another 5 ml of methylated spirit and grind for 3 minutes. Filter the solution.

The colours of the flower dyes in table 1 relate to a very suitable indicator, a cornflower extract.

The experiment can be used as an introduction to acids, only Evaluation - Question 3 in the report required a deeper understanding of the topic, it can also be used to introduce polyprotic acids.

Divide the students in groups, each group will conduct the experiment with only one indicator. The group's results should be exchanged among one another.

Disposal

After use, the solutions can be collected in the collecting tank for waste acids and bases for disposal.

Equipment

| Position No. | Material | Order No. | Quantity |
|---------------------|---|-----------|----------|
| 1 | Protecting glasses, clear glass | 39316-00 | 1 |
| 2 | Test tube rack, d = 22 mm, wood, for 6 test tubes | 37685-10 | 1 |
| 3 | Wash bottle, plastic, 250 ml | 33930-00 | 1 |
| 4 | Pipette with rubber bulb | 64701-00 | 6 |
| 5 | Laboratory pencil, waterproof | 38711-00 | 1 |
| 6 | Test tubes, d = 18 mm, l = 18 cm, glass, 100 pc | 37658-10 | 5 |
| | Potassium hydroxide, 500 g | 30103-50 | |
| | Litmus solution, 100 ml | 30127-10 | |
| | Sodium hydroxide, 500 g | 30157-50 | |
| | Hydrochloric acid, 37 %, 1000 ml | 30214-70 | |
| | Water, distilled, 5 l | 31246-81 | |
| | Acetic acid, 99-100%, 500 ml | 31301-50 | |
| | Methyl orange, 0.1 %, 250 ml | 31573-25 | |
| | Phenolphthalein, 0.5 % in ethanol, 100 ml | 31715-10 | |
| Additional material | | | |
| | Beetroot | | |
| | Flowers | | |
| | Natural dye (red cabbage) | | |



Safety information



Hazard and precautionary statements

Hydrochloric acid (5%)

| | |
|---------------------|--|
| H314: | Causes severe skin burns and eye damage. |
| H335: | May cause respiratory irritation. |
| H290: | May be corrosive to metals. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P301 + P330 + P331: | IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. |
| P309 +P310: | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |
| P305 + P351 + P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |

Acetic acid (10%)

| | |
|---------------------|--|
| H226: | Flammable liquid and vapour. |
| H314: | Causes severe skin burns and eye damage. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P301 + P330 + P331: | IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. |
| P307 + P310: | IF exposed: Immediately call a POISON CENTER or doctor/physician. |
| P305 + P351 + P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |

Sodium hydroxide (5%)

| | |
|---------------------|---|
| H314: | Causes severe skin burns and eye damage. |
| H290: | May be corrosive to metals. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P301 + P330 + P331: | IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. |
| P309 +P310: | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |
| P305 + P351 + P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses. Continue rinsing. |

Potassium hydroxide (5%)

| | |
|---------------------|---|
| H314: | Causes severe skin burns and eye damage. |
| H302: | Harmful if swallowed. |
| P280: | Wear protective gloves/protective clothing/eye protection/face protection. |
| P301 + P330 + P331: | IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Remove contact lenses. Continue rinsing. |
| P305 + P351 + P338: | IF IN EYES: Rinse cautiously with water for several minutes. |
| P309 +P310: | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |

Phenolphthalein

| | |
|-------|--|
| H226: | Flammable liquid and vapour. |
| P210: | Keep away from heat/sparks/open flames/hot surfaces. — No smoking. |

Hazards

- Acids and bases are irritating.
- Wear protecting glasses!

The effect of acids and bases on natural and synthetic indicators

(Item No.: P7510200)


Introduction

Application and task

How can you distinguish acids from bases?

Application

We encounter acids and bases in our everyday life as well as in the occasional chemistry lesson. They are present in the form of vinegar or in soaps. They are everywhere. In order to be able to handle acids and bases safely, it is important to identify those. One method to distinguish acids from bases is the use of indicators.



Experiment set-up

Task

Examine the different effects acids and bases have on natural as well as synthetic indicators.

Equipment



| Position No. | Material | Order No. | Quantity |
|---------------------|---|-----------|----------|
| 1 | Protecting glasses, clear glass | 39316-00 | 1 |
| 2 | Test tube rack, d = 22 mm, wood, for 6 test tubes | 37685-10 | 1 |
| 3 | Wash bottle, plastic, 250 ml | 33930-00 | 1 |
| 4 | Pipette with rubber bulb | 64701-00 | 6 |
| 5 | Laboratory pencil, waterproof | 38711-00 | 1 |
| 6 | Test tubes, d = 18 mm, l = 18 cm, glass, 100 pc | 37658-10 | (5) |
| | Potassium hydroxide, 500 g | 30103-50 | |
| | Litmus solution, 100 ml | 30127-10 | |
| | Sodium hydroxide, 500 g | 30157-50 | |
| | Hydrochloric acid, 37 %, 1000 ml | 30214-70 | |
| | Water, distilled, 5 l | 31246-81 | |
| | Acetic acid, 99-100%, 500 ml | 31301-50 | |
| | Methyl orange, 0.1 %, 250 ml | 31573-25 | |
| | Phenolphthalein, 0.5 % in ethanol, 100 ml | 31715-10 | |
| Additional material | | | |
| | Beetroot | | |
| | Flowers | | |
| | Natural dye (red cabbage) | | |

Set-up and procedure

Set-up

Hazards

- Acids and bases have a strong irritating effect!
- Wear protective glasses!



Set-up

Number the test tubes from 1 to 5 (Fig. 1).



Fig. 1

Put the test tubes next to each other in the test tube rack (Fig. 2) and label the pipettes (Fig. 3).



Fig. 2



Fig. 3

Procedure

Procedure

Use a pipette to fill test tube 1 with distilled water (fill level ca. 2 cm). Use the appropriate pipette (for each solution) to fill the test tubes with the same amount (fill level ca. 2 cm) as listed below:

- test tube 2: hydrochloric acid
- test tube 3: acetic acid
- test tube 4: sodium hydroxide solution
- test tube 5: potassium hydroxide solution

Use a new pipette to add a few drops of your chosen indicator to each test tubes. Note down your observations. Exchange your results with the other groups.



Fig. 4

Disposal

After use, the solutions can be collected in the collecting tank for waste acids and bases for disposal.

Report: The effect of acids and bases on natural and synthetic indicators

Result - Observations

Briefly describe your observations during the experiment (for example, the change of colour).

.....

.....

.....

.....

Result - Table 1

Complete the table below with your notes.

| Indicator | H ₂ O | | HCl | | C ₂ H ₄ O ₂ | | NaOH | | KOH | |
|-----------------|------------------|---|------------|---|--|---|-------------|---|-------------|---|
| Beetroot | wine red | 1 | brown red | 1 | brown red | 1 | brown green | 1 | brown green | 1 |
| Red cabbage | purple | 1 | red | 1 | red | 1 | blue green | 1 | blue green | 1 |
| Flower | dark purple | 1 | red | 1 | red | 1 | blue purple | 1 | blue purple | 1 |
| Litmus | dark purple | 1 | red | 1 | red | 1 | blue | 1 | blue | 1 |
| Methyl orange | yellow/orange | 1 | red | 1 | red | 1 | yellow | 1 | yellow | 1 |
| Phenolphthalein | colourless | 1 | colourless | 1 | colourless | 1 | pink | 1 | red purple | 1 |

Evaluation - Question 1

Make assumptions concerning the experiment and your observations.

.....

.....

.....

.....

Evaluation - Question 2

Which indicator is not suitable to identify acids?

.....

.....

.....

.....

Evaluation - Question 1

The indicator bromthymol blue is yellow when in acidic pH range and blue when in alkaline pH range. Which conclusion can you make, when the indicator shows a green colour?

.....

.....

.....

.....