
## Mechanical Equivalent of Heat Apparatus 04440.00

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### 1. APPLICATION

The mechanical equivalent of heat equipment provides a method of determining the mechanical equivalent of heat and of determining the specific thermal capacity of a solid body.

Mechanical energy from friction is fully converted into thermal energy. As, since 1977, both forms of energy are given in the same units of  $1\text{J} = 1\text{Nm}$ , the mechanical equivalent of heat = 1.

In the experiment a metal test body is rotated and heated by the frictional effect of a tensioned plastic band.

### 2. DESCRIPTION AND CONSTRUCTION

The equipment consists of a baseplate with a journal bearing in which the friction cylinder and the crank handle are screwed together with a holder for the friction band.

- Fix the friction cylinder and the crank handle to the journal bearing.

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- Hang a spring balance from the holder.
- Attach the friction band to the spring balance.
- Wrap the band 2.5 times around the cylinder so that the load on the balance is relieved if the crank handle is turned to the right.
- Attach a weight to the lower end of the friction band.
- Fill the hole in the cylinder with thermally conducting paste.

The temperature of the cylinder should be measured before and after a given number of rotations. The hole filled with thermally conducting paste is provided for this purpose. The temperature measurement is most accurate when the thermometer also has contact with the friction cylinder while it is being turned. The thermometer is held with an angled universal clamp with swivel (Fig. 1).

### IMPORTANT!

The setup is such that the hole in the cylinder and thermometer are precisely lined up. The baseplate must be clamped firmly with the screw clamps to the table top and must not move when the crank handle is turned. The thermometer may break if the apparatus moves and tilts.

### 3. PROCEDURE AND EVALUATION

The mass of the friction cylinder should be determined by weighing. A number of rotations, e.g. 100, are made with the crank handle as evenly as possible. The spring balance then indicates a constant force  $F_1$ . The work  $W$  due to friction on the cylinder is:

$$W = 2 \pi r n (F_2 - F_1)$$

$n$  = number of rotations

$r$  = cylinder radius

$F_1$  = force on the spring balance

$F_2$  = weight

The work is converted in thermal energy  $Q$ .

$$Q = c m (T_2 - T_1)$$

$m$  = mass of the friction cylinder

$c$  = specific thermal capacity of the cylinder

$T_1$  = cylinder temperature at start

$T_2$  = cylinder temperature at finish

For a known specific thermal capacity for the cylinder material (CuZn:  $c = 0.385 \text{ J/gK}$ ), the mechanical equivalent of heat  $W/Q$  can be found from the measurement.

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If it is assumed that  $W = Q$ , then the specific thermal capacity of solid bodies can be found with this apparatus. The experiment can be carried out with friction cylinders other than that supplied with the apparatus.

#### Note:

With an aluminium cylinder (Order No. 04441.03) the work due to friction must not be greater than 10 N (weight 1 kg) to avoid severe abrasion of the cylinder which might contaminate the friction band.

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#### 4. LIST OF EQUIPMENT

##### Included accessories:

|                                        |          |
|----------------------------------------|----------|
| 1 Friction cylinder CuZn, 0.64kg       | 04441.01 |
| 1 Plastic friction band                | 04441.04 |
| 1 Crank handle                         | 04441.05 |
| 1 Thermometer -10...+30°C, 0.1°C divs. | 05948.00 |
| 1 Thermally conducting paste, 50g      | 03747.00 |
| 2 Screw clamps                         | 02014.00 |

##### Accessories required:

|                       |          |
|-----------------------|----------|
| 1 Weight, 5kg         | 44096.81 |
| 1 Spring balance 100N | 03060.04 |

##### Holder for thermometer:

|                               |          |
|-------------------------------|----------|
| 1 Table clamp "Pass", small   | 03948.55 |
| 1 Universal clamp with swivel | 37716.00 |

##### For the determination of thermal capacity:

|                                  |          |
|----------------------------------|----------|
| 1 Friction cylinder CuZn, 1.28kg | 04441.02 |
| 1 Friction cylinder Al, 0.39kg   | 04441.03 |
| 1 Weight 1kg                     | 44096.70 |
| 1 Spring balance 10N             | 03060.03 |

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