



Neuro-simulator
Neuro-simulator, power supply

65963.00
65963.93

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Phone +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.de
Internet www.phywe.de


 The unit complies with the corresponding EC guidelines.



Fig. 1: Neuro-simulator 65963.00.

Fig. 2: Neuro-simulator, power supply 65963.93.

Operating instructions

1 SAFETY PRECAUTIONS



- Carefully read these operating instructions completely before operating this instrument. This is necessary to avoid damage to it, as well as for user-safety.
- Check that your mains supply voltage corresponds to that given on the type plate fixed to the instrument.
- Install the instrument so that the on/off switch and the mains connecting plug are easily accessible.
- Do not cover the ventilation slots.
- Take care that no liquids or objects enter in through the ventilation slots.
- Only use the instrument in dry rooms in which there is no risk of explosion.
- Do not start up this instrument should there be visible signs of damage to it or to the line cord.
- Only use the instrument for the purpose for which it was designed.

2 PURPOSE AND CHARACTERISTICS

The electronic neuron component, the "NEURO-SIMULATOR" enables the active simulation, i.e. defined by the user, of electrophysiological experiments from the cellular level up to the network level. This is achieved in that the intracellular potentials and the effects of the synapses can be derived as in the experiment on an "intracellular electrode". The same applies to the action potentials which can be derived through an "extracellular electrode" on the axon of the cell simulated by hardware. The cells are linked to one another and to a stimulation device with three buttons and an

optical sensor. Therefore, signals picked up via the axon can be applied to one or more synapses of the following "neurons" (neuronal network).

Each neuron component has a total of 9 synapses which can be connected. Here, axons (simulated by leads) of the sensors or other NEURO-SIMULATOR components can be terminated. An anatomical link is then possible. As in the nervous system, the effects of a signal that reaches a synapse of a following cell via an axon are solely determined by the synapse characteristics. The types of synapse are labelled by a colour code and can therefore be unambiguously identified.

An adjustable threshold defines which proportion of the intracellular depolarisation is passed on as action potentials via the efferent axon. The action potentials can be obtained at a signal tapping which symbolises an extra-cellular electrode. They can then be displayed on an oscilloscope or a personal computer with a suitable interface. In addition, an acoustic monitor, which can be deactivated, is included in each NEURO-SIMULATOR component so that the action potentials can be monitored audibly as is usual in electrophysiological research laboratories. The negative resting potential and the post-synaptic potentials can be obtained at the tapping point for the symbolised intracellular electrodes. The level of this potential can also be read from the light intensity of the light emitting diode at the tip of the electrode.

3 FUNCTIONAL AND OPERATING ELEMENTS

The neuro-simulator and the power supply are accommodated in plastic (ABS) housings. The cover plate of the housing has a carrying handle that can be swung upwards. The base plate has a similar handle, which can be swung out to enable the instrument to stand in an inclined




Fig. 3: Front view of the Neuro-simulators 65963.00.

positioned. Four rubber feet ensure slip-resistance and stability. The instrument can be stacked on other instruments having the same type of housing, whereby the rubber feet stand in the pan-shaped hollows of the instrument below for increased security against displacement. When units are stacked, the sloped position is only permissible for the topmost unit.

The unit is connected to the AC mains using the mains lead that is supplied with it. The lead is inserted into the equipment connection plug at the back of the unit.

The mains switch to start operating the unit is situated in the immediate vicinity of the equipment connection plug at the back of the unit. The neuro-simulator does not have its own power supply, but is instead supplied from the separate power supply. The front panels of the units are fitted with the following functional elements and controls (see Fig. 3 and 4):

NEURO-SIMULATOR

- 1 *Indicating lamp*
to indicate that the unit is switched on.
- 2 *Stimulator synapses*
(depolarising); 3 green 4 mm sockets for signal input.
- 3 *presynaptic synapses*
(quiescent inhibition); 3 brown 4 mm sockets for signal input.
- 4 *Hebb's synapse*
(variable depolarising) blue 4 mm socket for signal input.
- 5 *Inhibitory synapses*
(hyperpolarising).
2 red 4 mm sockets for signal input.
- 6 *Intracellular electrode (I)*
yellow 4 mm socket for the measurement of the resting potential and the post-synaptic potential.

- 7 *Illuminated display*
for displaying level of the post-synaptic potential.
- 8 *Control knob for threshold*
for setting the proportion of the intracellular depolarisation which is passed as the action potential via the efferent axon.
- 9 *Input/output*
for the supply voltage; DIN socket for the input of the $\pm 9V$ supply voltage produced in the power supply. The voltage can be looped through to another neuro-simulator using the second DIN socket.
- 10 *Ground*
white 4 mm socket as ground connection for intracellular (I) and extra-cellular (E) measurements.
- 11 *Efferent axon*
black 4 mm socket for obtaining the efferent (leading away) axon signal.
- 12 *Extracellular electrode (E)*
yellow 4 mm socket for the measurement of the action potential on the axon.
- 13 *Loudspeaker button*
for switching the loudspeaker on and off.
- 14 *Loudspeaker*
for the acoustic display of the action potentials.
- 15 *Reset button*
or resetting the Hebb's synapse to the basic state (fast "synaptic memory loss").




Fig. 4: Front view of the Neuro-simulator, power supply 65963.93.

NEURO-SIMULATOR, POWER SUPPLY

- 1 *Indicating lamp*
for indicating that the unit is switched on.
- 2 *Control knobs for stimulation intensity*
for setting the level of stimulation in stimulation channels 1, 2 and 3.
- 3 *Stimulation buttons*
A stimulation signal is present on the two corresponding sockets when the stimulation button is pressed.
- 4 *Stimulation outputs*
4 yellow sockets for measuring the stimulation intensity.
- 5 *Stimulation outputs*
4 black 4 mm sockets for obtaining the stimulation voltage for the synapses of a neuro-simulator.
- 6 *Ground*
2 white 4 mm sockets as ground connection for the measurement of stimulation intensity.
- 7 *Supply voltage output*
DIN socket for the output of the ± 9 V supply voltage for a maximum of 4 neuro-simulators.
- 8 *Offset voltage*
blue 4 mm socket for obtaining a voltage of -7 V. This is needed for the display of an intracellular potential on instruments not having any offset adjustment.
- 9 *Photoelectric sensor*
for picking up light stimuli (e.g. pocket torch). The output voltage present on the corresponding sockets is inversely proportional to the luminous intensity.

Changing the primary safety fuse

The fuse holder is in the upper part of the mains socket of the instrument, and so is only accessible when the connecting cord is not plugged in. Unplug the connecting cord, open the fuse holder using a screwdriver, take out the defect fuse and replace it with a new one (first check the specification of this against the data on the type plate), then fit the fuse holder back in the mains socket.

Should this fuse blow when the instrument is switched on, never replace it with a more resistant fuse! A defect is indicated and the instrument must be returned to the Phywe service department for repair.

4 NOTES ON OPERATION

This high-quality instrument fulfills all of the technical requirements that are compiled in current EC guidelines. The characteristics of this product qualify it for the CE mark.

This instrument is only to be put into operation under specialist supervision in a controlled electromagnetic environment in research, educational and training facilities (schools, universities, institutes and laboratories).

This means that in such an environment, no mobile phones etc. are to be used in the immediate vicinity. The individual connecting leads are each not to be longer than 2 m.

The instrument can be so influenced by electrostatic charges and other electromagnetic phenomena that it no longer functions within the given technical specifications. The following measures reduce or do away with disturbances:

Avoid fitted carpets; ensure potential equalization; carry out experiments on a conductive, earthed surface, use screened cables, do not operate high-frequency emitters (radios, mobile phones) in the immediate vicinity.

5 TECHNICAL SPECIFICATIONS (typical for 25°C)

Operating temperature range	5 - 40°C
Relative humidity	< 80%
NEURO-SIMULATOR	(Art.-No.: 65963.00)
Intracellular electrode (I)	
Resting potential	approx. -7 V (\pm -70 mV in the nerve cell)
Depolarisation (stimulation)	approx. 0 V (\pm 0 mV)
Hyperpolarisation (inhibition)	approx. -8.8 V (\pm -88 mV)
Extra-cellular electrode (E)	
Action potential (AP)	approx. 50 mV (\pm 500 μ V)
Duration of an AP	approx. 1 ms
Membrantime constant	approx. 50 ms (rise time) ca. 400 ms (decay time)
Supply voltage	\pm 9 V
Power consumption	2.5 W
Housing dimensions (mm)	230 \times 236 \times 236 (W, D, H)
Weight	approx. 2.8 kg

NEURO-SIMULATOR, POWER SUPPLY (Art.-No.: 65963.93)

Voltage output	\pm 9 V (for max. 4 Neuro-simulators)
Stimulator outputs	0...7 V
Offset output	-7 V
Mains supply	
Protection class	I
Connecting voltage (+6%/-10%)	see type plate
Mains frequency	50/60 Hz
Power consumption	14 VA
Mains fuse (5 mm \times 20 mm)	see type plate
Housing dimensions (mm)	230 \times 236 \times 168 (W, D, H)
Weight	approx. 2.6 kg

6 EXPERIMENT LITERATURE

Handbook Neuro-simulator	01191.01
Praktikumseinheit Nervenphysiologie	16703.41

7 ACCESSORIES

A set of connecting leads is included with the neuro-simulator:

Signal lead, l = 350 mm	11055.00
Connecting cord, 15 A, l = 6 cm, white	(2x) 168805
Connecting cord, 15 A, l = 15 cm, white	168806
Connecting cord, 15 A, l = 50 cm, white	07314.06
Connecting cord, 32 A, l = 50 cm, yellow	07361.02

8 LIST OF EQUIPMENT

Neuro-simulator	65963.00
Neuro-simulator, power supply	65963.93
Neuro-simulator, equipment set 4 + 1 (4 x 65963.00 + 1 x 65963.93)	65963.88

For displaying the experimental results:

Cobra3 BASIC-UNIT	12150.00
-------------------	----------

Power supply 12 VDC/2 A	12151.99
Data cable, plug/socket, 9 pole	14602.00
Software Cobra3 Universal recorder	14504.61
or:	
Oscilloscope 30 MHz, 2 channels	11459.95

PC, Windows®

9 GUARANTEE

We guarantee the instrument supplied by us for a period of 24 months within the EU, or for 12 months outside of the EU. This guarantee does not cover natural wear nor damage resulting from improper handling.

The manufacturer can only be held responsible for the function and technical safety characteristics of the instrument, when maintenance, repairs and changes to the instrument are only carried out by the manufacturer or by personnel who have been explicitly authorized by him to do so.

10 WASTE DISPOSAL

The packaging consists predominately of environmental compatible materials that can be passed on for disposal by the local recycling service.



Should you no longer require this product, do not dispose of it with the household refuse. Please return it to the address below for proper waste disposal.

PHYWE Systeme GmbH & Co. KG
Abteilung Kundendienst
Robert-Bosch-Breite 10
D-37079 Göttingen

Phone +49 (0) 551 604-274
Fax +49 (0) 551 604-246