Printed: 13.04.2017 12:01:42 | P1001900

Buoyancy and floating (Item No.: P1001900)

Curricular Relevance Experiment: Subtopic: Area of Expertise: **Education Level: Topic:** Mechanik der Auftrieb und Physik Klasse 7-10 Mechanik Flüssigkeiten Schwimmen Difficulty **Preparation Time Execution Time Recommended Group Size** 22222 00000 \odot 2 Students 10 Minutes 10 Minutes Easy **Additional Requirements: Experiment Variations: Keywords:**

Task and equipment

Information for teachers

Additional Information

The students should determine whether or not various objects float in water. They should then investigate whether the volume of an object or its shape affect its buoyancy.

Remarks

Since the plasticine mass should not be specified, the number of little balls that the ship can carry is also unspecified. The important thing is that the students realise that when a certain mass is exceeded, the carrying capacity of the ship becomes insufficient and the ship sinks.



PHYWE

advanced



Buoyancy and floating (Item No.: P1001900)

Task and equipment

Task

When does an object float?

You will immerse various solid objects with different densities in water and determine whether or not they float. Then you will form a hollow object from a plastic material (plasticine) and determine if it floats.





advanced

Equipment



Position No.	Material	Order No.	Quantity
1	Aluminium column	03903-00	1
1	Wood column	05938-00	1
2	Glass beaker DURAN®, short form, 600 ml	36015-00	1
2	Plasticine, 10 sticks	03935-03	(1)
2	Rubber stopper 26/32, 1 hole 7 mm	39258-01	1
2	Tubing connect.,T-shape,ID 8-9 mm	47519-03	1
3	Rubber ball,diam.15 mm	03921-00	1

Printed: 13.04.2017 12:01:42 | P1001900



Set-up and procedure

- Fill the beaker three-quarter full of water and place the wood and aluminium columns (Fig. 1, 2), the rubber ball, the piece of plasticine, the connecting piece (plastic) and the rubber stopper one after another in the water.
- Observe whether the objects sink or float; record your observations under "Result Observations 1" in the report.



- Divide the piece of plasticine into two pieces of about equal size. Form a ball from one half and a thin sheet from the other half from which you should shape a boat (Fig. 3).
- Place both pieces in water and observe whether they sink or float (Fig. 4). Tick the appropriate checkboxes in "Result Table 1", part (a), in the report.
- Now form a ball from the boat and a boat from the ball. Place both objects in the water again and repeat the experiment, ticking the appropriate checkboxes in "Result Table 1", part (b), in the report.





• Divide the plasticine ball into about 10 little balls which are as nearly as possible the same size. Place the boat in the water and load it bit by bit with the little balls. Record your observations under "Result - Observations 2" in the report.



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107



Report: Buoyancy and floating

Result - Observations 1

Sort the investigated material in floating and non-floating material.

Result - Table 1

Note down your observations in the table:

Object		Floats? (yes/no)	
(a)	Ball	no ¹	
	Boat	yes 1	
(b)	Ball	no 1	
	Boat	yes 1	



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 13.04.2017 12:01:42 | P1001900



Result - Observations 2

Note down your observations of the boat loaded with little balls:

Evaluation - Question 1

Which of the six objects floats?



Robert-Bosch-Breite 10 D - 37079 Göttingen

.....

Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 13.04.2017 12:01:42 | P1001900



Evaluation - Question 2

Is the buoyancy of the individual objects dependent on the material they are made of?

Evaluation - Question 3

Is there a correlation with the object's densities? Arrange the objects according to whether their density is larger or smaller than that of water ($\rho = 1 \text{ g/cm}^3$) and compare these results with your recordings under "Result - Observations 1". What can you conclude?



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 13.04.2017 12:01:42 | P1001900



Evaluation - Question 4

Is the buoyancy of the objects dependent on their volume?

Evaluation - Question 5

Draw a floating unloaded and a floating loaded boat in water. What is the difference?

.....



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 13.04.2017 12:01:42 | P1001900



Evaluation - Question 6

Why does the plasticine ball sink, while the boat made of the same material and with the same mass does not sink?

.....

Evaluation - Question 7

What is the carrying capacity of the boat?



Robert-Bosch-Breite 10 D - 37079 Göttingen