

Forced oscillation and resonance (Item No.: P1003100)

Curricular Relevance



Difficulty



Intermediate

Preparation Time



10 Minutes

Execution Time



10 Minutes

Recommended Group Size



2 Students

Additional Requirements:

Experiment Variations:

Keywords:

Task and equipment

Information for teachers

Additional Information

In the Experiment "Helical spring pendulum" the students have determined the natural frequency of a spring pendulum and established this frequency dependence on pendulum mass and spring constant. In this experiment they should force a spring pendulum to oscillate and from observations of the forced oscillation amplitude make conclusions concerning the correlation between excitation frequency and natural frequency.

Furthermore, the natural frequency of a spring pendulum should be measured and then compared with the (approximate) natural frequency determined in the second part of the experiment.


Forced oscillation and resonance (Item No.: P1003100)

Task and equipment


Task

How can the oscillations of a spring pendulum be induced?

Start a spring pendulum oscillating with your hand and observe the effect.
Measure the oscillation frequency with which the spring pendulum oscillates normally.



Equipment



Position No.	Material	Order No.	Quantity
1	Support base, variable	02001-00	1
2	Support rod, stainless steel, l = 600 mm, d = 10 mm	02037-00	1
3	Boss head	02043-00	1
4	Holding pin	03949-00	1
5	Helical spring, 3 N/m	02220-00	1
6	Set of precision weights, 1g-50g	44017-00	1
7	Stop watch 4	03078-00	1
Additional material			
	Scissors		

Set-up and procedure

Set-up

First screw the splitted support rod together (Fig. 1). Set up a stand with the support base (Fig. 2), put the support rod in the support base and tight it with the screw (Fig. 3).



Fig. 1




Fig. 2




Fig. 3

Fix the bosshead to the support rod. Fix the holding pin in the bosshead and hang the helical spring in it (Fig. 4). Attach a 50 g mass piece from the weight set to the helical spring (Fig. 5).



Fig. 4



Fig. 5

Procedure

- Deflect the spring pendulum downwards and let it oscillate at its natural oscillation frequency (Fig. 6).
- Start the stop watch when the pendulum is at its lower reversal point and measure the time required for 10 complete oscillations.
- Repeat this measurement twice and record the times in Table 1 in the report.



Fig. 6

- Take hold of the upper end of the helical spring (Fig. 7).
- Move your hand up and down with the spring pendulum very slowly (low excitation frequency). Observe the movement of the spring pendulum and record your observations in the report.
- Move your hand more rapidly than before (intermediate excitation frequency) and again observe the spring pendulum.
- Move your hand more rapidly (high excitation frequency, i.e. larger than the natural oscillation frequency) and again observe the spring pendulum.




Fig. 7

In order to disassemble the support base you should press the yellow buttons (Fig. 8).



Fig. 8

Report: Forced oscillation and resonance

Results - Table 1

Record all the measured values in Table 1.

Calculate the average value for 10 oscillations and from it the time for one oscillation, i.e. the oscillation period T .

Use the oscillation period to calculate the pendulum's natural oscillation frequency of $f_0 = 1/T$.

Measurement N.	t_{10} in s	Average t_{10} in s	T in s	f_0 in Hz
1	1	1	1	1
2	1	1 ± 0	1 ± 0	1 ± 0
3	1	1 ± 0	1 ± 0	1 ± 0

Results - Observations 1

Write down your observations at low excitation frequency.

.....

.....

.....

.....

Results - Observations 2

Write down your observations at intermediate excitation frequency:

.....

.....

.....

.....

Results - Observations 3

Write down your observations at high excitation frequency:

.....

.....

.....

.....

Evaluation - Question 1

How is the amplitude affected by

- a) low excitation frequency?
- b) intermediate excitation frequency?
- c) high excitation frequency?

.....

.....

.....

.....

Evaluation - Question 2

How can you explain the observed behavior?

.....

.....

.....

.....

Evaluation - Question 3

One uses the term "resonance" to describe the case in which the excitation frequency and the natural frequency of an object capable of oscillation agree. How is the resonance of a spring pendulum noticeable?

.....

.....

.....

.....

Evaluation - Question 4

What determines the oscillation frequency of a spring pendulum?

.....

.....

.....

.....

Evaluation - Question 5

Who is the "exciter" in this case?

.....

.....

.....

.....