Linear expansion of metals (Item No.: P1042900)



### Task and equipment

#### Information for teachers

### **Additional Information**

Metal tubes are heated with steam. The linear expansion of three different metals is determined with the aid of a rotating-shaft pointer. From this the metals' coefficient of linear expansion is determined (Supplementary problem). The use of a rolling-shaft pointer is advantageous because the resulting indication of expansion is frictionless. The functional principle of a rotating-shaft pointer can be easily demonstrated in a simple experiment and studied by the students (see Question 4). The students should do this before they calculate the change in length  $\Delta I$  in "Evaluation".

#### Remarks

- 1. The setup must not be touched while the tube is being heated; otherwise, the sensitive pointer will be displaced.
- 2. No condensed water must accumulate in front of or in the metal tube; otherwise the tube will not be heated to a uniform 100 °C.
- 3. The pointer's deflection must be marked as an extension of the pointer (see Fig. 22).

### Notes to the angle determination

The angle  $\varphi$ , which is formed by the pointer's deflection s, can best be determined using the scale in the figure. This scale only applies for a pointer length (to the table surface) of l = 10.5 cm. The angle can also be calculated with the formula

 $\tan \varphi = s/l.$ 

For this calculation the distance  $\Delta l/2$  is neglected with reference to *s*.

The students should test the functional principle of the rotating-shaft pointer, e.g. by rolling it through an angle of 360° with a ruler (see Fig. 3), and by doing so determine the following: Movement of the ruler  $\Delta I = 2.5$  cm.

Movement of the rotating-shaft pointer  $\Delta l = 1.25$  cm. In general:

 $\Delta l/2 = 2\pi \times r \times \varphi/360^\circ.$ 



Robert-Bosch-B D - 37079 Götti Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 **HYWE** 

### **Teacher's/Lecturer's Sheet**

Printed: 24/07/2019 17:39:59 | P1042900







Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107



# Linear expansion of metals (Item No.: P1042900)

### Task and equipment

#### Task

### How does the expansion of metals change during heating?

Measure the expansion of steel, brass and aluminium while heating them from room temperature to 100 °C.





## Equipment



| Position No.           | Material  | Order No. | Quantity |
|------------------------|---|-----------|----------|
| 1                      | Support base, variable                              | 02001-00  | 1        |
| 2                      | Support rod, stainless steel, I = 250 mm, d = 10 mm | 02031-00  | 1        |
| 2                      | Support rod, stainless steel, I = 600 mm, d = 10 mm |           | 2        |
| 3                      | Boss head   |           | 3        |
| 4                      | Ring with boss head, i. d. $= 10$ cm                | 37701-01  | 1        |
| 4                      | Universal clamp                                     | 37715-00  | 1        |
| 5                      | Wire gauze with ceramic, 160 x 160 mm               | 33287-01  | 1        |
| 6                      | Beaker, low form, plastic, 100 ml                   | 36011-01  | 1        |
| 7                      | Erlenmeyer flask 100 ml, wide-neck SB 29            | 36428-00  | 1        |
| 8                      | Glass tube, straight, l=80 mm, 10/pkg.              | 36701-65  | (1)      |
| 8                      | Students thermometer,-10+110°C, l = 180 mm          | 38005-02  | 1        |
| 9                      | Rubber stopper 26/32, 1 hole 7 mm                   | 39258-01  | 1        |
| 10                     | Measuring tape, I = 2 m                             | 09936-00  | 1        |
| 10                     | Silicone tubing i.d. 7mm                            | 39296-00  | 1        |
| 11                     | Collar for linear expansion                         | 04231-55  | 1        |
| 12                     | Brass tube  | 04234-01  | 1        |
| 12                     | Aluminium tube                                      | 04234-03  | 1        |
| 12                     | Iron tube   | 04234-02  | 1        |
| 13                     | Rotating shaft with pointer                         | 04236-01  | 1        |
|                        | Butane burner, Labogaz 206 type                     | 32178-00  | 1        |
|                        | Butane cartridge C206, without valve                | 47535-01  | 1        |
|                        | Glycerol, 250 ml                                    | 30084-25  | 1        |
|                        | Boiling beads, 200 g                                | 36937-20  | 1        |
| Additional<br>material |   |           |          |
|                        | Matches   |           |          |
|                        | Felt-tip pen  |           | 1        |
|                        | Sheet of paper (approx, 100 mm x 50 mm)             |           | 1        |

Robert-Bosch-Breite 10 D - 37079 Göttingen

advanced

### Set-up and procedure

### Set-up

#### Warning!

- 1. Always insert the thermometer or glass tubes into the rubber stopper using glycerol.
- 2. During the heating of the water the support ring and the wire gauze get extremely hot!
- 3. The metal tube is hot, steam escapes at its end!
- 4. The setup must not be touched during the course of the experiment; otherwise, the sensitive pointer will be displaced.

#### Setup

• Set up the support stand according to the following pictures.





Printed: 24/07/2019 17:39:59 | P1042900





Fig. 4

Fig. 5



Fig. 6

Fig. 7



• Place the collar for linear expansion in the notch (spout side) on the aluminium tube.



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 24/07/2019 17:39:59 | P1042900





• Attach the tube on the other side of the location of the notch in the bosshead.



- The metal tube should be slightly slanted so that condensed water vapour (steam) can run out.
- For this reason attach the bosshead on the spout side so that it touches the collar.



Fig. 12

• Place the rotating shaft with pointer between the bosshead and the collar.

Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107





• Select the height of the metal tube so that the pointer tip is as close to the table surface as possible.



• Place the beaker under the tube's spout.







• Fill the Erlenmeyer flask half-full with water and add two beads.



Fig. 16

Fig. 17

• Insert the glass tube into the stopper and seal the flask.

Printed: 24/07/2019 17:39:59 | P1042900





Fig. 18

Fig. 19



Fig. 20

• Position the tubing in such a way that no condensed water accumulates in front of the tube.



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107

Printed: 24/07/2019 17:39:59 | P1042900



### Procedure

- Under the pointer's tip tape a piece of paper (approx. 5 cm x 10 cm) to the table.
- Set the pointer perpendicular and mark the initial position of the pointer's tip.



Fig. 21

- Measure the pointer's length from its pivot point to the table top. This distance / should be 10.5 cm. If it is not, correct your experimental setup.
- Measure the room temperature  $\vartheta_0$  and record it in the report.
- Bring the water to a boil.
- Wait until condensed water and steam escape from the tube.
- Watch the pointer until it stops moving.
- Mark the pointer's new position (distance *s* in Fig. 2: use your ruler to determine where the point would now touch the paper).
- Turn off the burner, wait until the tube has cooled somewhat and then repeat the experiment with other tubes (remove the hot tube with a cloth).



advanced

**PHYWE** 

## **Report: Linear expansion of metals**

#### **Result - Observations 1**

| 1. | Length of the metal tube:             | $l_0 = 50 \text{ cm}$                 |    |
|----|---------------------------------------|---------------------------------------|----|
| 2. | Temperature of the heated tube:       | $\vartheta_1 = 100 \ ^\circ \text{C}$ |    |
| 3. | Radius of the pointer's rotating shaf | t: $r = 0.2 \text{ cm}$               |    |
| 4. | Length of the pointer (to table top): | l = 10,5  cm                          |    |
| 5. | Room temperature:                     | $\vartheta_0 =$                       | ٥С |

#### Result - Table 1

- 1. For each of the three metal tubes measure the pointer's deflection s on the sheet of paper and record them in the table.
- 2. Determine the angle  $\varphi$  (in degrees) for each pointer deflection *s* (in cm) with the aid of the scale in the figure 23 and record it in the table.
- 3. When the metal tube expanded by the length  $\Delta l$ , the rotating pointer moved a distance  $\Delta l/2$  further. The arc of angle is therefore as large as the distance  $\Delta l/2$ :

#### $\Delta l/2 = 2\pi \cdot r \cdot \phi/360^{\circ}$

Calculate the linear expansion  $\Delta l$  for the different metal tubes and record it in the table.

| Material  | <i>s</i> in cm | φ in Degrees | ∆/ in cm |
|-----------|----------------|--------------|----------|
| Aluminium | 1              | 1            | 1        |
|           | ±0             | ±0           | ±0       |
| Brass     | 1              | 1            | 1        |
|           | ±0             | ±0           | ±0       |
| Iron      | 1              | 1            | 1        |
|           | ±0             | ±0           | ±0       |

Printed: 24/07/2019 17:39:59 | P1042900



#### **Evaluation - Question 1**

Compare the metals in terms of their expansion, i.e. arrange them in order of decreasing expansion.

#### **Evaluation - Question 2**

Compare the thermal expansion of metals to that of liquids or gases.



Robert-Bosch-Breite 10 D - 37079 Göttingen

Printed: 24/07/2019 17:39:59 | P1042900



#### **Evaluation - Additional Task 1**

The linear expansion  $\Delta l$  is a function of the temperature difference  $(\vartheta_1 - \vartheta_0)$  and of the total length  $l_0$ . Could this be shown with your experimental setup?

#### **Evaluation - Additional Task 2**

The linear expansion of rods and tubes is characterised by the coefficient of linear expansion  $\alpha$ . It is expressed by the equation:

#### $\Delta l = \alpha \, l_0 \, (\vartheta_1 - \vartheta_0).$

Calculate the coefficient  $\boldsymbol{\alpha}$  and record it in the table.

|           | Literature value          |                           |
|-----------|---------------------------|---------------------------|
| Material  | α in 10 <sup>-6</sup> /°C | α in 10 <sup>-6</sup> /°C |
| Aluminium | 23,7                      | 23                        |
| Brass     | 18,3                      | 17                        |
| Iron      | 10                        | 11                        |



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107