


Determination of density by measurement of the buoyancy (Item No.: P1297300)

Curricular Relevance



Difficulty



Intermediate

Preparation Time



10 Minutes

Execution Time



10 Minutes

Recommended Group Size



1 Student

Additional Requirements:

- Spirit, ca. 400 ml

Experiment Variations:

Keywords:

Principle and equipment

Principle

It is to be shown, how the density of solid and liquid substances can be determined from a measurement of the buoyancy.

Equipment

| Position No. | Material | Order No. | Quantity |
|----------------------|-----------------------------------|-----------|----------|
| 1 | Demo Physics board with stand | 02150-00 | 1 |
| 2 | Torsion dynamometer | 03069-03 | 1 |
| 3 | Sinker, aluminum | 03903-01 | 1 |
| 4 | Support plate on fixing magnet | 02155-00 | 1 |
| 5 | Glycerol, 250 ml | 30084-25 | 2 |
| 6 | Glass beaker DURAN®, tall, 600 ml | 36006-00 | 1 |
| Additional material: | | | |
| | Spirit, ca. 400 ml | | |

Set-up and procedure

1. Experiment

- Fasten the torsion dynamometer at the upper edge of the demo-board.
- Hang the sinker on, measure the weight force of it in air ($F_{G,L}$) and note this value.
- Position the shelf underneath the sinker at the lower edge of the board, then place the glass beaker containing about 400 ml of water on it (Fig. 1).
- Lower the dynamometer with sinker, until the sinker is completely immersed in the water.
- Measure the weight force $F_{G,W}$ with which the sinker now pulls on the dynamometer, and note this value.




Fig. 1

2. Experiment

- Use the set-up as in Experiment 1, but pour alcohol into the glass beaker instead of water.
- Lower the dynamometer, with the sinker hung on it, until the sinker is completely immersed in the alcohol.
- Measure the weight force $F_{G,Sp}$ with which the sinker pulls on the dynamometer and note this value.
- Replace the alcohol in the beaker with glycerol and carry out the same procedure to measure and note the force $F_{G,Gl}$.

Observation and evaluation

Observation

Evaluation