

# **Changeover switches and alternate switches**

(Item No.: P1380500)

#### **Curricular Relevance**



## **Principle and equipment**

## Principle

This demonstration shows how switching between two electrical devices is possible, and how alternate switching functions.







Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 Printed: 25/10/2017 11:30:30 | P1380500



## Equipment

Position No.	Material	Order No.	Quantity
1	PHYWE power supply, universal DC: 018 V, 05 A / AC: 2/4/6/8/10/12/15 V, 5 A	13500-93	1
2	Demo Physics board with stand	02150-00	1
3	Switch, change-over, module DB	09402-02	2
4	Socket for incandescent lamp E10 ,module DB	09404-00	2
5	Junction, module DB	09401-10	2
6	Electr.symbols f.demo-board,12pcs	02154-03	1
7	Connector, straight, module DB	09401-01	2
8	Connector, angled, module DB	09401-02	4
9	Connector, T-shaped, module DB	09401-03	1
10	Filament lamps 12V/0.1A, E10, 10 pieces	07505-03	1
11	Connecting cord, 32 A, 1000 mm, red	07363-01	1
12	Connecting cord, 32 A, 1000 mm, blue	07363-04	1



Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 info@phywe.de www.phywe.com



## Set-up and procedure

#### **1st. Experiment**

- Connect up the circuit as shown in Fig. 1, the power supply is set to 0 V and switched off.
- Switch on the power supply and adjust the voltage to the 12 V rated voltage of the lamp; observe the lamp (1).
- Repeatedly operate the changeover switch and observe the behaviour of the lamp (2).
- Switch off the power supply.



#### **2nd. Experiment**

- Connect up the circuit as shown in Fig. 2.
- Switch on the power supply and repeatedly operate the changeover switches in any way you want, observe the lamp (3).
- Set the power supply to 0 V and switch it off.









Robert-Bosch-Breite 10 D - 37079 Göttingen Tel: +49 551 604 - 0 Fax: +49 551 604 - 107 info@phywe.de www.phywe.com



### **Observation and evaluation**

#### Observation

- 1. Only one lamp lights up when the circuit is closed.
- 2. The lamps light up alternately when the changeover switch is switched.
- 3. The lamp can be switched on and off with each of the two changeover switches.

#### **Evaluation**

A changeover switch enables you to switch over in a circuit from one electrical appliance to another one. With two changeover switches and three connections (one of these to the source of current), you can close or break a circuit from two different positions any way you wish. This is called alternate switching.

#### Remarks

Alternate switching is mostly used in the lighting of staircases and long corridors. The changeover switches which are used for alternate switching are called alternate switches.

